

MgCl₂ (25mM)

#GP012.1250 (1.25 ml) | GP012.6250 (6.25 ml)
 (FOR RESEARCH ONLY)



Product: Magnesium ion (Mg²⁺) is an essential cofactor for DNA polymerase activity, as it enables the incorporation of dNTPs during DNA polymerization. Moreover, Mg²⁺ facilitates the binding of primers to DNA templates by stabilizing negative charges, thereby influencing annealing and melting temperatures.

High Mg²⁺ concentrations increase polymerase activity, however, at the cost of fidelity. dsDNA is stabilized by elevated Mg²⁺ concentrations thereby preventing denaturation and increasing the melting temperature. For the same reason, primer-dimer formation is increased with higher Mg²⁺ concentrations as well as annealing temperatures.

Bioavailability of Mg²⁺ in a PCR reaction is dependent on many parameters such as template DNA concentration and purity, primer concentration and purity, dNTP concentration, presence of chelators (e.g. EDTA) or the presence, type and concentration of fluorescent dyes. Thus, it is recommended to optimize the Mg²⁺ concentration for each new PCR reaction or whenever changing these parameters.

Applications: PCR

Quantity: #GP012.1250 contains 1.5 ml of 25mM MgCl₂ in PCR-grade water.
 #GP012.6250 contains 5 vials of 1.5 ml of 25mM MgCl₂ in PCR-grade water.

Appearance: Clear transparent liquid

Storage: Store at 2-8°C for up to 6 months or at -20°C for up to 2 years.

Optimization:

The normal range of Mg²⁺ concentration used in PCR is 0.5mM to 5mM. Magnesium ions can be added either in the form of MgCl₂ or Mg₂SO₄. A normal starting concentration is 1.5mM MgCl₂. It is common practice that DNA polymerases are supplied with reaction buffers that already contain MgCl₂ in such a concentration that the final concentration will be 1.5mM. In the table hereunder the amount of 25mM MgCl₂ to be added to a reaction mixture without supplemented MgCl₂ as well to a reaction mixture already containing MgCl₂ (1.5mM final concentration) is given:

Final [MgCl ₂] in a:	0.5mM	1.0mM	1.5mM	2.0mM	2.5mM	3.0mM	3.5mM	4.0mM	4.5mM	5.0mM
50µl reaction mixture w/o supplemented MgCl ₂	1.0	2.0	3.0	4.0	5.0	6.0	7.0	8.0	9.0	10.0
25µl reaction mixture w/o supplemented MgCl ₂	0.5	1.0	1.5	2.0	2.5	3.0	3.5	4.0	4.5	5.0
20µl reaction mixture w/o supplemented MgCl ₂	0.4	0.8	1.2	1.6	2.0	2.4	2.8	3.2	3.6	4.0
50µl reaction mixture with supplemented MgCl ₂	x	x	0.0	1.0	2.0	3.0	4.0	5.0	6.0	7.0
25µl reaction mixture with supplemented MgCl ₂	x	x	0.0	0.5	1.0	1.5	2.0	2.5	3.0	3.5
20µl reaction mixture with supplemented MgCl ₂	x	x	0.0	0.	0.8	1.2	1.6	2.0	2.4	2.8